



SPECIFICATIONS

SUGAR SPHERES SANAQ[®] (USP40-NF35 / Ph. EUR)

Type	250-355	355-500	500-600	600-710	710-850	850-1000	1000-1180	Test Method	
Definition	Sugar Spheres are white nearly spherical pellets with a slight sweet taste.								
Identification								A B C	Ph. EUR Ph. EUR Ph. EUR
Particle Size								In House	
	0 % > 425 ≤ 10 % > 400 ≥ 90 % > 300 ≤ 10 % < 300	0 % > 600 10 % > 500 90 % > 355 10 % < 355	0 % > 710 10 % > 600 90 % > 500 10 % < 500	0 % > 850 10 % > 710 90 % > 600 10 % < 600	0 % > 1000 10 % > 850 90 % > 710 10 % < 710	0 % > 1180 10 % > 1000 90 % > 850 10 % < 850	0 % > 1400 10 % > 1180 90 % > 1000 10 % < 1000		
Friability								n.m.t. 5.0 %	Ph. EUR / USP
Sucrose Content								62.5 - 91.5 %	Ph. EUR / USP
Loss on drying								n.m.t 4.0 %	Ph. EUR / USP
Heavy Metals								< 5 ppm	Ph. EUR / USP
Sulphate Ash /Residue on ignition								0.20%	Ph. EUR / USP
Microbiological Parameters									
TAMC								10 ² CFU/g	Ph. EUR / USP
TYMC								10 ² CFU/g	Ph. EUR / USP
E. coli, Pseudomonas aeruginosa, Staphylococcus Aureus								negative in 1 g sample	Ph. EUR / USP
Salmonella species								negative in 10 g sample	Ph. EUR / USP

Contact

Pharmatrans Sanaq AG

Gewerbestrasse 18
CH-4123 Allschwil - Basel
Switzerland
Tel.: +41 61 225 90 00
Email: info@pharmatrans-sanaq.com
Web: www.pharmatrans-sanaq.com